MFT – BODYTEAMWORK

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INSTRUCTION MANUAL

Intended Usage

MFT training devices can be used in athletic and fitness training to improve ones coordination and proprioception and in physical coordination, therapy for injury prevention and rehabilitation of ones ankles, knees, hips and spine.

Requirements

Constant supervision and instruction from a doctor, certified trainer or therapist is required when using MFT training devices for physical therapy.

Age

from 4 to 100 years, depending on training level

Weight

max. 120 kg body weight (80 kg Fun Disc, 100 kg Fun Disc Digital)

Skills

no special skills required

Safety Instructions

One must execute all workouts on a compact, non slippery, dry and level surface. The device must be clean and dry during all training. Wear sneakers or train barefoot. Do not wear socks, stockings, etc.. One can reduce the risk of injuries when falling by keeping the workout area clear of debris and if needed, place a soft surface around the training area. Children must only exercise under supervision of an adult.

Additional Instruction

MFT devices must not be used for jumping, but only for smooth and easy movements.

Regulations/Maintenance

Before usage check the device for damages and make sure that all adjustable pieces are fixed.

www.mft-bodyteamwork.com
office@mft-bodyteamwork.com

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Activation of MFT Balance Sensor

- ▶ Activate Bluetooth on the end device (tablet, smartphone, computer).
- ▶ Download and start the app/software and move/tap the disc briefly. If several devices are detected, select the appropriate device.

Calibration

- ▶ The disc must be on the floor when starting the app. No one is allowed to be on the disc!
- ▶ The disc is positioned correctly when the MFT logo faces towards you.
- ▶ Do not change the position of the disc during training, otherwise it will have to be recalibrated.

MFT Bodyteamwork App
Test/Training/Games

MFT Volare App
Gaming

MFT Training & Workouts

Batteriewechsel / Change Battery

MFT DIGITAL LINE

DIN EN 50419

Nach/according to: 2012/19/EU

WEEE-Reg.Nr.: DE65255934 (ElektroG)

Reg.Nr.: 46489055 (BattG)

UN-Nr. 3091