
MFT BALANCE SENSOR SIT BALL

Batteriewechsel / Change Battery / Replacement de la batterie / Sostituzione della batteria
MFT Balance Sensor: Varta CR 2032 LiMn 3V
Fehleranalyse & Videoanleitungen
Troubleshooting & More:
www.mft-bodyteamwork/support

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Software-download / mehr Infos:

www.mft-bodyteamwork.com/support

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Important information:

Software-download / more information::

www.mft-bodyteamwork.com/en/support

App download: AppStore, GooglePlay, Amazon Appstore

Commissioning

• Remove the battery fuse strip before first use (photo 1)

Activation of MFT Balance Sensor with Bluetooth

• Activate Bluetooth on the device (tablet, smartphone, computer)

• Move the ball slightly – tap – start app/software.

If multiple devices are detected, select the appropriate device.

Calibration & seating position

• The ball must be on the ground when starting the app

• Sit on the middle of the ball with both feet on the ground

• The balance sensor is correctly positioned when the MFT logo is

pointing forward - between the legs and towards the terminal (see photo 2)

• During training, the location of the ball should not be changed, otherwise it

must be recalibrated.

System requirements:

Bluetooth from 4.0 LE, Apple from iOS 6, Mac OSX from 10.7, Android from 4.3,

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
MFT Trainings Devices Instruction Manual

Intended Usage: MFT training devices can be used in athletic and fitness training to improve ones coordination and proprioception and in Physical Therapy for injury prevention and rehabilitation of ones ankles, knees, hips and spine.

Requirements: Constant supervision and instruction from a doctor, certified trainer or therapist is required when using MFT training devices for Physical Therapy.
Age: from 4 to 100 years, according to the level of fitness.
Weight: max. 120 kg body weight (approx. 270 lbs)
Skills: no special skills necessary

Safety instructions: One must execute all workouts on a compact, non slippery, dry and level surface must be clean and dry during all training. The product cannot be combined with other products. This product must not be used as a toy. The device must be clean for use and must not have a slippery surface. Changes to the device are not permitted. Ensure that there is no free fall space and, if necessary, place fall-resistant underlays on top. For safety reasons, we generally recommend the use of a ball bowl. Children may only train under the supervision of an adult person.

Additional information: The device is not suitable for jumping or standing, but only for slow, sensitive movements while sitting.

Regulations/Maintenance: Before usage check the device for damages.

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