Balance  
Trainer  
2.0

# CHALLENGE DISC 2.0

// GEBRAUCHSANWEISUNG // INSTRUCTION MANUAL  
// MODE D'EMPLOI // ISTRUZIONI D'USO  
// GUIDA ALL'ALLENAMENTO

[www.mft-bodyteamwork.com](http://www.mft-bodyteamwork.com)

  
**MFT**<sup>®</sup>  
Bodyteamwork

# CHALLENGE DISC 2.0

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**MFT Bodyteamwork App |**

Test/Training/ Games

**MFT Volare App | Gaming**

**Software-download:**

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**Mehr Infos:**

[www.mft-bodyteamwork.com/support](http://www.mft-bodyteamwork.com/support)

## Important information:

**MFT Bodyteamwork App |**

Test/Training/ Games

**MFT Volare App | Gaming**

**Software-download:**

[www.mft-bodyteamwork.com/support](http://www.mft-bodyteamwork.com/support)

**App download:**

AppStore, GooglePlay, Amazon Appstore

### Activation of MFT Balance Sensor with Bluetooth

- Activate Bluetooth on the device (tablet, smartphone, computer)
- Move the disc slightly – tap – start app/software. If multiple devices are detected, select the appropriate device.

### Calibration

- The disc must be on the ground when starting the app. You must not stand on the disc at this stage!
- The disc/balance sensor is correctly positioned when the MFT logo is pointing forward – between the legs and towards the terminal
- During the training, the location of the ball should not be changed, otherwise it must be recalibrated.

### Option: MFT App on TV

- Tablet/smartphone – screen synchronization via Apple TV Box, Google Chrome, Amazon Fire TV, etc.

**More information:**

[www.mft-bodyteamwork.com/support](http://www.mft-bodyteamwork.com/support)

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[www.mft-bodyteamwork/support](http://www.mft-bodyteamwork/support)

**Remarque  
importante relative:**

**MFT Bodyteamwork App |**

**Test/Training/ Games**

**MFT Volare App | Gaming**

**Software-download:**

[www.mft-bodyteamwork.com/support](http://www.mft-bodyteamwork.com/support)

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UN-Nr. 3091

<https://www.mft-bodyteamwork.com/CE/>

## MFT Trainings Devices Instruction Manual

**Intended Usage:** MFT training devices can be used in athletic and fitness training to improve ones coordination and proprioception and in Physical Therapy for injury prevention and rehabilitation of ones ankles, knees, hips and spine.

**Requirements:** Constant supervision and instruction from a doctor, certified trainer or therapist is required when using MFT training devices for Physical Therapy.

**Age:** from 4 to 100 years, according to he level of fitness.

**Weight:** max. 120 kg body weight (aprox: 270 lbs)

**Skills:** no skills necessary

**Safety Instructions:** One must execute all workouts on a compact, non slippery, dry and level surface must be clean and dry during all training. One can train with sneakers or barefoot. Do not wear socks, stockings or dress shoes. If necessary, one can use a chair, table, bar, or spotting person for assistance during fitness training. One can reduce the risk of injuries when falling by keeping the workout area clear of debris and if needed, place a soft surface around the training area. Children must only exercise under supervision of an adult.

**Additional Instruction:** MFT devices must not be used for jumping, but only for smooth and easy movements.

**Regulations/Maintenance:** Before usage check the device for damages and make sure that all adjustable pieces are fixed.

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